

Drug Checking Best Practices

Fentanyl and xylazine test strips

(03/2026)

Fentanyl test strips

Use only with stimulants, like cocaine, methamphetamine, MDMA. **Substances where fentanyl is not expected nor desired.**

Xylazine test strips

Use only with opioids, like fentanyl, heroin, and prescription pills. **False positives are common when checking stimulants for xylazine.**

Test strip results

Read the test in bright light. Faint lines are REAL lines. Two red lines means negative. One red line means positive.

CAUTION: The accuracy of drug checking technologies is not guaranteed. The safest practice is to avoid using drugs altogether.

Steps to use test strips



dhs.wi.gov/drug-checking

Negative

Positive

