

Rapid Testing Temperature Log

Instructions

Document storage temperatures of both test kits and the controls for all business days. Additional instructions and details can be found in the [HIV Counseling, Testing and Referral Services Protocol](#) under the "Ensuring Proper Temperatures for Tests and Controls" section.

Thermometer location: _____

Month/Year (XX/XXXX): _____

Acceptable temperature range*: _____

Day	Initials	High temperature	Low temperature	Corrective action taken when temperature is out of range
1				
2				
3				
4				
5				
6				
7				
8				
9				
10				
11				
12				
13				
14				
15				
16				
17				
18				
19				
20				
21				
22				
23				
24				
25				
26				
27				
28				
29				
30				
31				

*The acceptable range for **Determine test kit** storage is 2° to 30° C or 36° to 86° F and the acceptable range for control storage is 2° to 8° C or 36° to 46° F. The acceptable range for the **Rapid Syphilis Health Check Test** is 4° to 30° C or 39° to 86° F and the acceptable range for syphilis control storage is 2° to 8° C or 36° to 46° F.

Reviewed by (print name): _____

Date reviewed: _____