



Pediatric Equipment Checklist

Mannitol or 3% saline
Intubation, code and neurologic medications
Catheter-over-the-needle device 22 gauge
Catheter-over-the-needle device 24 gauge
Pediatric intraosseous needles or device
Intravenous solutions: normal saline
Intravenous solutions: dextrose 5% normal saline
Infant C-collars
Child C-collars
Cuffed endotracheal tube: 3.5 millimeters
Cuffed endotracheal tube: 4.5 millimeters
Cuffed endotracheal tube: 5.5 millimeters
Cuffed endotracheal tube: 6.5 millimeters
Laryngoscope: Straight: 1
Laryngoscope: Straight: 2
Laryngoscope: Curved: 2
Infant nasopharyngeal airways
Child nasopharyngeal airways
Oropharyngeal airway, size 0
Oropharyngeal airway, size 1
Oropharyngeal airway, size 2
Oropharyngeal airway, size 3
Oropharyngeal airway, size 4
Pediatric stylets for endotracheal tubes
Infant suction catheters
Child suction catheters
Bag-mask device, self-inflating: infant: 450 milliliters
Masks to fit bag-mask device adaptor for infants and children
Clear oxygen masks: partial non-breather infant
Clear oxygen masks: partial non-breather child
Infant nasal cannulas
Child nasal cannulas

	Nasogastric tubes: Infant: 8 French size
	Nasogastric tubes: Child: 10 French size
	Laryngeal mask airway: size 1.5
	Laryngeal mask airway: size 2
	Laryngeal mask airway: size 2.5
	Laryngeal mask airway: size 3
	Chest tubes: Infant: 10 or 12 French size
	Chest tubes: Child: one in the 16-24 French size range